

IKA Werke GmbH & Co. KG

CalWin

Overview


Content

1	Starting CalWin.....	4
2	User interface – basics.....	5
2.1	Hide and expand areas.....	5
2.1.1	Vertical splitters.....	5
2.1.2	Horizontal splitter.....	8
2.2	Search filter.....	9
2.2.1	Selecting or defining a search criteria.....	10
2.3	Main menu buttons.....	11
3	CalWin first steps.....	12
3.1.1	XML.....	12
3.1.2	SQL.....	12
3.2	Splash Screen.....	13
3.3	Initial Main Screen.....	14
3.4	Setup a calorimeter (C1, C200, C6000 iso, C6000 gs).....	15
3.4.1	Setup a calorimeter (C5000, C2000).....	17
3.4.2	Deactivate a calorimeter.....	17
3.4.3	Delete a calorimeter.....	18
3.4.4	Change calorimeter settings.....	18
3.5	Setup the calorimeter vessels.....	19
3.5.1	Automatic vessel import.....	19
3.5.2	Manual vessel input.....	21
3.5.3	Filtering Vessels.....	24
3.5.4	Delete or change a vessel.....	24
3.6	Import Measurements.....	25
4	Prepare a new measurement.....	28
4.1	Run a measurement.....	30
4.2	Cancel running measurement.....	33
4.3	Measurement finished.....	34
4.4	Evaluating the measurement.....	34
4.5	Make modifications to an existing measurement.....	36
4.6	States of a measurement.....	38
5	Settings.....	39
5.1	User Management.....	39
5.1.1	Create a new user.....	40
5.2	External Energies.....	43
5.3	C-Values.....	48
5.3.1	Comparing C-Values.....	49


5.4	System Settings	51
5.4.1	Grid Columns 1 – 3.....	52
5.4.2	Database Settings	53
5.5	Groupings	54
6	Detail tabs.....	57
6.1	Details	57
6.2	Graph	58
6.3	Analysis.....	58
6.4	Notes.....	59
6.5	Printing/Export Details	60

1 Starting CalWin

After you have successfully installed CalWin on your PC, you can start it from the Windows Start Menu



Or from the desktop





2 User interface – basics


2.1 Hide and expand areas

The areas of the main screen can be adjusted by the user.

2.1.1 Vertical splitters


You can press on the left splitter to hide the search area of the main screen










If you want to reshove the search area, press on the splitter to expand.

The width other vertical splitter can be adjusted by grapping the splitter (left mouse click) and move to the desired width or by clicking on the splitter to hide.






2.1.2 Horizontal splitter



2.2 Search filter

To restrict the shown measurements in the “Selected Measurements” area of the main list it is possible to define one or more search criteria for a desired search.





The following search example should show the possibility of a search:

- All measurements performed **yesterday** with **calorimeter xxx** from **user xxx** that are in the **group xxx** and have the substring **calib** in their sample name.

2.2.1 Selecting or defining a search criteria


Press on a button to expand the search tree and select the criteria by radio button. The selected value is shown in the button. To start the search, press on the “Update Measurement List” button.



2.3 Main menu buttons



1. **New Measurement:** Opens the dialog for setting up a new measurement.
2. **Settings:** Opens the dialog for the system settings
3. **Export:** Performs an EXCEL export of the currently selected measurements with all the columns that are visible. To change the visible columns to get more information in the EXCEL file, read the chapter 5.4.1 about setting up the visible columns.
4. **Print:** Opens the printer setup dialog to print the currently selected measurements with their visible columns.
5. **Help:** Activates the help system and shows additional tooltip information.




Press the “Help” button again to deactivate this option.

INFORMATION: If a button is hidden by its tooltip, click on it and the tooltip closes.

The help system is currently not implemented for the whole application! At the moment it only includes the items on the main screen!

3 CalWin first steps

When opening CalWin for the first time the following screen appears to select the data storage process.



Here you can select between:

- Storing locally in a XML file. This method is recommended for most users!
- Using a SQL Database

3.1.1 XML

Stores the data (measurements, vessels, c-values, users,...) in a structured file in the user's documents folder.

3.1.2 SQL

Uses the Microsoft SQL Server to store the data. Therefore the Microsoft SQL Server 2008 R2 must be installed. The advantage of the SQL server is that it can be installed locally on your PC or on a network server. For a detailed installation description please read the following document:

SQL Server 2008 for CalWin.pdf

The setup process can be found in this document:

Setup Database.pdf

You can find these documents included on the DVD.


3.2 Splash Screen

The splash screen displays the software version and the loading process. As soon as all modules are fully loaded, CalWin's main screen is shown.




3.3 Initial Main Screen

The first step to start working with CalWin is to setup a new calorimeter and its vessel(s). Afterwards you can create a new measurement and let it run.



3.4 Setup a calorimeter (C1, C200, C6000 iso, C6000 gs)

Click on „Settings“ in the main menu and select the Calorimeter button.



1. Start the calorimeter and select its main menu view.
2. Connect the calorimeter to a free PC COM-port.
3. Click on „New Calorimeter“
4. Select the COM-port, the calorimeter is connected to.



5. Click on „Calorimeter Information“ to receive the serial number and calorimeter type automatically from the device.

Calorimeter Information

Calorimeter Serial Number: 00.0000000

Calorimeter Name:

Calorimeter Type: C1

Measurement Method: Isoperibol 22°C

COM-Port: COM3

In Use

Add Calorimeter
Delete Calorimeter
Change Calorimeter
Reset Calorimeter

Calorimeter Information

6. Enter a calorimeter name and click on "Add Calorimeter".

Calorimeter Information

Calorimeter Serial Number: 00.0000000

Calorimeter Name: C1 No 1

Calorimeter Type: C1

Measurement Method: Isoperibol 22°C

COM-Port: COM3

In Use

Add Calorimeter
Delete Calorimeter
Change Calorimeter
Reset Calorimeter

Calorimeter Information

7. Now the calorimeter appears in the overview list below.

Calorimeter Information

Calorimeter Serial Number: 00.0000000

Calorimeter Name: C1 No 1

Calorimeter Type: C1

Measurement Method: Isoperibol 22°C


COM-Port: COM3

In Use

New Calorimeter
Delete Calorimeter
Change Calorimeter
Reset Calorimeter
Import Vessels
Import Measurements
Calorimeter Information

Name	Calorimeter Type	Measurement Mode	COM Port	SampleRack	Calorimeter ID	EqualCalori
C1 No 1	C1	Isoperibol 22°C	COM3		00.0000000	

8. Now the calorimeter appears in the main view of CalWin



3.4.1 Setup a calorimeter (C5000, C2000)

Here it is not possible to receive the calorimeter information automatically. All inputs have to be done manually:


1. **Calorimeter Serial Number:** This can be found on a label at the backside of the device.
2. **Calorimeter Name:** The name of the device (e.g. "C5000 No.1")
3. **Calorimeter Type:** The type of the device (C5000, C2000)
4. **Measurement Method:** The calorimeter working mode.
5. **The COM-Port:** The port your calorimeter is connected to.

If all settings are done, press on the "Add Calorimeter"-button to add the calorimeter.


3.4.2 Deactivate a calorimeter

When a calorimeter is currently out of use, you can deactivate it by unchecking the "In Use" box located under "Settings" > "Calorimeter".

To access the „In Use“ checkbox, click on "Change Calorimeter" and then check or uncheck the checkbox. Finally click "Save Changes" to adopt the new value.




After unchecking „In Use“, the calorimeter disappears from the main screen and the COM port is free to be used with a different calorimeter or device.
 Use this option, if you just want to temporarily work without this calorimeter.




3.4.3 Delete a calorimeter

If you want to remove a calorimeter from your CalWin software, you can do so by clicking „Delete Calorimeter“, while accessing the calorimeter settings.



3.4.4 Change calorimeter settings

If you want to change the name, the measurement method or the COM port of a calorimeter, you can do so by clicking „Change Calorimeter“. The parameters that are subject to change get activated. After finishing your changes, click “Save Changes”.




3.5 Setup the calorimeter vessels

3.5.1 Automatic vessel import

It is possible to import the vessels including their c-values directly from the device. This works for the C1 and C6000 calorimeter. You find the option to import the vessel under „Settings“ > „Calorimeters“. As soon as connection to a calorimeter is established on the main form, the option “Import Vessels” is activated. The c-values are imported from the calorimeter device.

This is the recommended procedure for the C6000iso and C6000gs Calorimeter, because the vessel’s serial number must also be the RFID-tag and can be read automatically by the device.



The imported vessels can now be found under „Vessels“.

Settings

User Management
External Energies
Calorimeter
Vessels
C-Values
System Settings
Groupings

Vessel Informations:

Vesse-ID: 1
Name: V1
Optical ID:
Calorimeter: C1 No 1 (Isoperibol 22°C | COM5)
Working Mode: Isoperibol 22°C
C-value [J/K]: 4111
Ignitions: 0
Note: Vessel imported from device C1 (C1 No 1)

New Vessel
Delete Vessel
Change Vessel
Reset Vessel
Delete vessel in all modes

Vessel Name	Vessel-ID/RFID	Optical ID	Measurement Mode	C-Value	Ignition Counter
V1	1		Isoperibol 22°C	4111	0
V2	2		Isoperibol 22°C	0	0
V3	3		Isoperibol 30°C	0	0
V4	4		Isoperibol 30°C	0	0

Search Calorimeter: All Calorimeters
Search Working Mode: All Measurement Modes
Search


Print | OK

3.5.2 Manual vessel input

The other option to enter the vessel information is to use the Vessel settings.

The screenshot shows the 'Settings' application window with the 'Vessel Informations' section. The 'New Vessel' button is highlighted with a red circle. The form includes input fields for Vessel-ID, Name, Optical ID, Calorimeter, Working Mode, C-value [J/K], Ignitions, and Note. Below the form is a table with columns: Vessel Name, Vessel-ID/RFID, Optical ID, Measurement Mode, C-Value, and Ignition Counter. At the bottom, there are search filters for Calorimeter and Working Mode, and buttons for Print and OK.

1. Press the "New Vessel" button.
2. Enter the Vessel-ID/RFID. This number can be found on the C6000gs/iso display.
This value must be unique!
Use the RFID number for the C6000gs/iso vessels and the engraved vessel serial number for other devices. The automatic vessel detection of the C6000 device only works with a correct RFID number.
3. Enter a vessel name
4. The optical ID can be set for the C5000/C2000 if the vessel detection is activated
5. Select a calorimeter from the calorimeter combo box.
6. Select the working mode of the calorimeter
7. Select a c-value (if already known) or leave initial value.
8. Enter the number of ignitions that have already been performed with this vessel.
9. **Optional:** It's possible to enter an additional note to the vessel, e.g., "in use since 01.01.2013")
10. Click on „Add Vessel“ to complete the adding process.



11. Now the vessel appears in the overview list.

Settings

- User Management
- External Energies
- Calorimeter
- Vessels**
- C-Values
- System Settings
- Groupings

Vessel Informations:

Vessel-ID: 11

Name: C1 Vessel No 1

Optical ID:

Calorimeter: C1 No 1 (Isoperibol 22°C | COM5)

Working Mode: Isoperibol 22°C

C-value [J/K]: 4115

Ignitions: 12

Note:

New Vessel

Delete Vessel

Change Vessel

Reset Vessel

Delete vessel in all modes

Vessel Name	Vessel-ID/RFID	Optical ID	Measurement Mode	C-Value	Ignition Counter
C1 Vessel No 1	11		Isoperibol 22°C	4115	12

Search Calorimeter: All Calorimeters

Search Working Mode: All Measurement Modes

Search

Print | OK

3.5.3 Filtering Vessels

In case you have a lot of vessels added and stored, it is possible to search for:

- All vessels of a specific calorimeter
- All vessels in a specific measurement mode
- All vessels of a calorimeter in a certain measuring mode

1. Select the calorimeter (optional)
2. Select the working mode (optional)
3. Press the “Search” button to start the search and update the vessel list.

3.5.4 Delete or change a vessel

By pressing the “Delete Vessel” or “Change Vessel” button, the user can delete or change an existing vessel.

- **Delete:** The vessel is no longer in use.
- **Change:** possible to change the name, optical ID (used for C5000, C2000), c-value, ignition counter or the note of a vessel.

The screenshot shows the 'Settings' application window. On the left is a sidebar with navigation options: User Management, External Energies, Calorimeter, **Vessels**, C-Values, System Settings, and Groupings. The main area is titled 'Vessel Informations:' and contains several input fields: Vessel-ID (11), Name (C1 Vessel No 1), Optical ID (empty), Calorimeter (C1 No 1 (Isoperibol 22°C | COM5)), Working Mode (Isoperibol 22°C), C-value [J/K] (4115), Ignitions (12), and Note (empty). To the right of these fields are buttons: New Vessel, Delete Vessel, Change Vessel, Reset Vessel, and Delete vessel in all modes. Below the input fields is a table with the following data:


Vessel Name	Vessel-ID/RFID	Optical ID	Measurement Mode	C-Value	Ignition Counter
C1 Vessel No 1	11		Isoperibol 22°C	4115	12

At the bottom of the window, there is a search section with two dropdown menus: 'Search Calorimeter:' (set to 'C1 No 1 (Isoperibol 22°C | COM5)') and 'Search Working Mode:' (set to 'Isoperibol 22°C'). A 'Search' button is to the right of these dropdowns. This search section is circled in red. At the very bottom are 'Print' and 'OK' buttons.

3.6 Import Measurements

It is possible to import measurements from the calorimeter device C6000gs/iso and C1 to CalWin. To perform this task, ensure that the calorimeter is connected and listed in the main menu (calorimeter must be in **WAIT** state).

Additional buttons (Import Vessels and Import Measurements) appear in the calorimeter settings menu.



1. Press the "Import Measurements" button
2. Wait until the user message disappears.

IKA® CalWin

Settings

Calorimeter Information

Calorimeter Serial Number: 00.0000000

Calorimeter Name: C1 No 1

Calorimeter Type: C1


Measurement Method: Isoperibol 22°C

COM-Port: COM5

In Use

Buttons: New Calorimeter, Delete Calorimeter, Change Calorimeter, Reset Calorimeter, Import Vessels, Import Measurements


Archive measurement transmission

 Please wait while measurement archive is imported from device C1 No 1!


Name	Calorimeter Type	Measurement Mode	COM Port	Calorimeter ID	In Use	Sample Rack	Decou
C1 No 1	C1	Isoperibol 22°C	COM5	00.0000000	<input checked="" type="checkbox"/>	<input type="checkbox"/>	

Print OK

3. The imported measurements and calibrations are now visible on the main screen.



4. All imported measurements are marked in the note field as "Imported from device"



5. It is not possible to see the graph of an imported measurement/calibration! The temperature values are not available if the measurement was not controlled by CalWin.

4 Prepare a new measurement

The screenshot displays the IKA CalWin 2012 software interface. The title bar at the top reads 'IKA® CalWin 2012'. The main menu bar includes 'New Measurement' (highlighted with a red circle), 'Settings', 'Export', 'Print', and 'Help'. Below the menu bar, the 'Selected Measurements:' section contains a table with columns: Calorimeter Name, Probe Name, Vessel Name, Weight [g], Result C-Value, Result (Ho), Vessel C-Value, and User Name. The table is currently empty. On the left side, there is a 'Period' dropdown menu set to 'All Measurements', and a 'Select Date:' section with 'from:' and 'to:' dropdowns both set to '12.02.2013'. Below this is a list of filters for 'Calorimeters', 'User', 'Groups', and 'Names'. At the bottom left, there is an 'Update Measurement List' button. The bottom right section shows a 'Device' window for 'C1 No 1' with a 'Waiting...' progress bar and a play button. To the right of the device window is a 'Details' panel with tabs for 'Details', 'Graph', 'Analysis', and 'Note'. The 'Details' tab is active, showing fields for 'Probe Name', 'Weight [g]', 'Result', 'User Name', 'Calorimeter Name', 'Vessel Name', and 'Start Date', all with 'init' values. Below these are 'External Energies' with 'Ext. Energy 1 Calculated Value', 'Ext. Energy 2 Calculated Value', 'Ext. Energy 3 Calculated Value', and 'Ignition Energy', all with 'init' values.

New measurement

Measurement

Measurement name: Calibration

User: Simulation

Calorimeter:

Vessel:

Weight [g]:

Group:

	Type:	Ho [J/g]:	Weight [g]:	QExt [J]:
<input checked="" type="checkbox"/> Combustion Aid 1	<input type="text" value="IKA Cotton Thread"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="50"/>
<input type="checkbox"/> Combustion Aid 2	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> Combustion Aid 3	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

Sum External Energies [J]:

Information:

1. Measurement name

Can be filled in by the user or left blank. If left blank, the system automatically fills in the current date and time as measurement name.

2. User

Users can be managed under „Settings“ > „User Management“. Users already created can be selected from the drop down list. It is also possible to enter a user name directly.

3. Calorimeter

Choose an available calorimeter from the drop down menu.

4. Vessel

The assigned vessels of the selected calorimeter are listed in this drop down menu.

5. Weight

The weight can either be entered manually or be transmitted by a connected balance.

6. Group

Groups can be managed in the settings under “Groupings”. It is possible to assign the measurement to a specific group to make it easier to find later on. The group can also be assigned if the measurement is finished.

7. Combustion Aids

Combustion Aids can be managed under „Settings“ > „External Energies“. These values can be selected from the drop down menu.

8. Calibration

Check this option to perform a calibration.

9. Simulation

Set this option to simulate the measurement/calibration. The user has to enter a ΔT value and the ignition energy.


4.1 Run a measurement

When all entries are complete, click on „OK“ to add the measurement to the main list.


The screenshot shows a software dialog box titled "New measurement". It contains several input fields and checkboxes. The "Measurement name" is "Test C1", "User" is "ttt (test user)", "Calorimeter" is "C1 No 1 (Isoperibol 22°C | COM5)", "Vessel" is "C1 Vessel No 1 (ID: 01.54678111 | Isop)", and "Weight [g]" is "1,0111". There are checkboxes for "Calibration" and "Simulation", both of which are unchecked. A "Copy Measurement" button is visible. Below these fields, there are three rows for "Combustion Aid" (1, 2, and 3). "Combustion Aid 1" is checked and has a type of "IKA Cotton Thread". The "Ho [J/g]" and "Weight [g]" fields for "Combustion Aid 1" are both "0", and the "QExt [J]" field is "50". The "Sum External Energies [J]" field is also "50". At the bottom, there are "Clear", "OK", and "Cancel" buttons. The "OK" button is circled in red.


Type:	Ho [J/g]:	Weight [g]:	QExt [J]:
<input checked="" type="checkbox"/> Combustion Aid 1 IKA Cotton Thread	0	0	50
<input type="checkbox"/> Combustion Aid 2			
<input type="checkbox"/> Combustion Aid 3			

Sum External Energies [J]: 50




If the calorimeter is in the Waiting state, it is possible to assign the measurement to the calorimeter by clicking on the assign button of the measurement or by holding down the right mouse button and dragging the measurement to the device panel before releasing the button.





When the measurement is assigned, the calorimeter changes its status from **Waiting** to **Prepared**. A prepared measurement can also be removed by clicking on the remove button.

If the calorimeter doesn't start automatically like the C1 calorimeter, you have to press the "start" button on the device or in CalWin



4.2 Cancel running measurement

The screenshot displays the IKA CalWin software interface. At the top, there is a menu bar with options: New Measurement, Settings, Export, Print, and Help. The main window is divided into several sections:


- Left Panel:** A sidebar with a 'Period' dropdown set to 'Today'. Below it are radio buttons for 'Today', 'Yesterday', 'Last Week', 'Last Month', 'Last Year', 'All Measurements', and 'Select Date:'. There are also date input fields for 'from: 21.02.2013' and 'to: 21.02.2013'. At the bottom of the sidebar are expandable sections for 'Calorimeters', 'User', 'Groups', and 'Names', along with an 'Update Measurement List' button.
- Selected Measurements:** A table with columns: Calorimeter Name, Probe Name, Vessel Name, Weight [g], Result C-Value, Result (Ho), Vessel C-Value, and User Name. The first row shows: C1 No 1, Test C1, C1 Vessel No 1, 1.0101, 0, 0, 5186. A red arrow points to the 'C1 No 1' entry with the text 'Measurement is running'.
- Device View:** A window titled 'C1 No 1' showing 'Test C1' with a progress bar and the text 'Closing (Fill O₂)'. A red arrow points to a 'Cancel' button in the top right corner of this window with the text 'Press here to cancel a running measurement'.
- Details Panel:** A panel on the right showing measurement details: Probe Name (Test C1), Weight [g] (1.0101), Result (Ho) (0), User Name (--), Calorimeter Name (C1 No 1), Vessel Name (C1 Vessel No 1), and Start Date (2/21/2013 3). Below this is an 'External Energies' section with values for Ext. Energy 1, 2, and 3 Calculated Values, and Ignition Energy (100).

There are two different possibilities when canceling a measurement:

- If the measurement is cancelled before the ignition takes place, it is possible to reassign and restart it.
- If cancelled after the ignition, the measurement is marked as cancelled and cannot be reassigned.


4.3 Measurement finished

When the measurement is finished, the result is displayed and the calorimeter changes back to the **Waiting** state.



4.4 Evaluating the measurement

When a measurement is finished, it can be accessed by double clicking to enter additional values needed for the selected report.



Selected Measurement: Test C1

Measurement Base | **Brennwert Korrektur** | Heizwert Korrektur | Ergebnisse

Measurement

Measurement name: Calibration

User: Simulation

Calorimeter: ΔT :

Vessel: IgnEnergy:

Weight:

Group:

User Defined Reference Combustion Value

Demo Mode:

Typ:	Ho:	Weight:	QExt:
<input checked="" type="checkbox"/> Combustion Aid 1: <input type="text" value="KA Zündfaden"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="50"/>
<input type="checkbox"/> Combustion Aid 2: <input type="text"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>
<input type="checkbox"/> Combustion Aid 3: <input type="text"/>	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>

Sum External Energies:

Result (Hov): 30115

Information:

Now the measurement menu displays additional tabs needed for the report.

Selected Measurement: Test C1

Measurement Base | **Brennwert Korrektur** | Heizwert Korrektur | **Ergebnisse**

Version: DIN 51900 1989

$Hov = (CVal * DT - (QN + QS + QCI + QF + QZ)) / mp$

CVal =

DT =

$QN = 0,97 * NO_2$ $NO_2 =$

$QS = 94,62 * S * mp$

Aus Schwefelgehalt S [%] $S =$ $mp =$

Aus entstandener Schwefelsäure

$QCI = 20,8 * CI * mp$ $CI =$ $mp =$

$QF = 32,4 * F * mp$ $F =$ $mp =$

$QZ = QIGN + QEXT1 + QEXT2 + QEXT3$

$QIGN =$

$QEXT1 =$

$QEXT2 =$

$QEXT3 =$

Ergebnisse:

$QN = 0$ $QS = 0$ $QCI = 0$ $QF = 0$ $QZ = 59$

Hov = 30115

When you are done click on „Evaluation Done“ to lock the data into place and save them from further manipulation. With „Print“ you are able to print out the form. Now the measurement is marked as „evaluated“ in the main list.

Selected Measurements:							
<input type="checkbox"/>			Calorimeter Name	Probe Name	Vessel Name	Weight [g]	Resu
<input type="checkbox"/>			C1 No 1	Test C1	C1 Vessel No 1	1.0101	0

4.5 Make modifications to an existing measurement

If you want to do some changes to an existing measurement, you first have to make a copy of this measurement as it is not allowed to edit the original data.

To do so, access to the original data and click on “Copy Measurement”.

The copy can now be found in the main list and is marked as a simulation.

Selected Measurement: Test C1 [X]

Measurement Base | Brennwert Korrektur | Heizwert Korrektur | Ergebnisse

Measurement

Measurement name: Calibration

User: Simulation

Calorimeter: ΔT :

Vessel: IgnEnergy:

Weight:

Group:

User Defined Reference Combustion Value

Demo Mode

Typ:

Combustion Aid 1 Ho: Weight: QExt:

Combustion Aid 2 Ho: Weight: QExt:

Combustion Aid 3 Ho: Weight: QExt:

Sum External Energies:

Result (Hov): 30115


Information:

The screenshot shows the IKA CalWin software interface. At the top, there is a menu bar with options: New Measurement, Settings, Export, Print, and Help. The main window is titled "Selected Measurements:" and contains a table with the following columns: Calorimeter Name, Probe Name, Vessel Name, Weight [g], Result C-Value, Result (Ho), Vessel C-Value, and User Name. The table has two rows. The first row is highlighted in blue and contains: C1 No 1, Test C1_0, C1 Vessel No 1, 1.0101, 0, 30115.47, 5186. The second row is also highlighted in blue and contains: C1 No 1, Test C1_1, C1 Vessel No 1, 1.0101, 0, 30115.47, 5186. A red circle is drawn around the second row, and a red arrow points to the simulation icon (a dollar sign with a slash) in the first column of that row. Below the table, there is a section for "Device" showing "C1 No 1" and a "Waiting..." status. To the right, there is a "Details" panel with tabs for "Details", "Graph", "Analysis", and "Note". The "Details" tab is active, showing the following information: Probe Name: Test C1_1, Weight [g]: 1.0101, Result Simulation Measurement: 30115.47, User Name: --, Calorimeter Name: C1 No 1, Vessel Name: C1 Vessel No 1, Start Date: 2/21/2013 3:38:03. Under the "External Energies" section, the values are: Ext. Energy 1 Calculated Value: 50, Ext. Energy 2 Calculated Value: --, Ext. Energy 3 Calculated Value: --, and Ignition Energy: 9.

Calorimeter Name	Probe Name	Vessel Name	Weight [g]	Result C-Value	Result (Ho)	Vessel C-Value	User Name
C1 No 1	Test C1_0	C1 Vessel No 1	1.0101	0	30115.47	5186	
C1 No 1	Test C1_1	C1 Vessel No 1	1.0101	0	30115.47	5186	

Copy of an original measurement marked as simulation


If you open the simulated measurement, it is possible to change the values of ΔT , ignition energy and you can change a calibration into a measurement and vice versa.







It is also possible to delete a simulated measurement. A regular measurement cannot be deleted.

4.6 States of a measurement

In the main list the measurements, calibrations and simulations are marked as such.



The meaning of the icons:

-  calibration
-  simulated calibration
-  measurement
-  simulated measurement

5 Settings

5.1 User Management

Settings

User Management

External Energies

Calorimeter

Vessels

C-Values

System Settings

Groupings

User Informations

Initials:

Name:

E-Mail:

Administrator

Password:

Confirm Password:

New User

Delete User

Change User

Cancel


Available Users

User Initials	User Name	Administrator	eMail
---------------	-----------	---------------	-------


Print | OK

5.1.1 Create a new user

CalWin can create and save various users. Measurements can be assigned to specific users.



1. Press "New User" button



2. Enter user initials, full name and/or e-mail address.
INFO: e-mail addresses are currently not in use. Later on it will be possible to send a lost password to an user!
3. Activate the administrator checkbox for the password entry fields.

INFO: administrator password is currently not in use! Later on, the Settings menu can be locked by password, so that a „normal“ user can only perform measurements without access to the settings.

The screenshot shows the 'Settings' dialog box with the 'User Management' section selected. The 'User Informations' form is filled with the following data:

Initials:	xyz
Name:	abc xyz
E-Mail:	abc@xyz.com
<input checked="" type="checkbox"/> Administrator	
Password:	****
Confirm Password:	****

The 'Add User' button is circled in red. Below the form is an 'Available Users' table:

User Initials	User Name	Administrator	eMail
---------------	-----------	---------------	-------

At the bottom of the dialog are 'Print' and 'OK' buttons.

4. When you are finished click on „Add User“ to complete the task.

The screenshot shows the 'Settings' dialog box with the 'User Informations' form empty. The 'Administrator' checkbox is unchecked. The 'New User' button is circled in red. Below the form is an 'Available Users' table:

User Initials	User Name	Administrator	eMail
xyz	abc xyz	<input checked="" type="checkbox"/>	abc@xyz.com

At the bottom of the dialog are 'Print' and 'OK' buttons.

5. From now on, this user will be available for selection when setting up a new measurement.

New measurement

Messung

Measurement name: Calibration

User: Simulation

Calorimeter:

Vessel:

Weight:

Group:

User Defined Reference Combustion Value

Demo Mode

	Typ:	Ho:	Weight:	QExt:
<input type="checkbox"/> Combustion Aid 1	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> Combustion Aid 2	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> Combustion Aid 3	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

Sum External Energies:

Information:

Copy Measurement

Clear OK Cancel

5.2 External Energies

It's possible to manage 3 different external energies. These energies can later be selected when setting up a new measurement. They can either have fixed energy or have a gross calorific value, which needs an additional weight indication.

The screenshot shows the 'Settings' window with the 'External Energies' section. The interface includes a sidebar with navigation options: User Management, External Energies, Calorimeter, Vessels, C-Values, System Settings, and Groupings. The main area displays the configuration for an external energy. At the top, three tabs are labeled 'External Energy 1', 'External Energy 2', and 'External Energy 3', with red circles and numbers 1, 2, and 3 above them respectively. The 'External Energy 1' tab is selected. The form contains the following fields and buttons:

- Name: IKA Cotton Thread
- Gross Calorific Value [J/g]: (empty)
- External Energy (fix) [J]: 50
- Note: Auto generated entry
- Buttons: New External Energy, Delete External Energy, Change External Energy, Reset External Energy

Below the form is a table with the following data:

External Energy Name	Gross Heat Value	External Energy
IKA Cotton Thread	0	50

At the bottom of the window are 'Print' and 'OK' buttons.

1. Select an external energy tab 1 – 3 to insert item.
2. Click on “New External Energy”

The screenshot shows the 'Settings' window with the 'External Energies' section active. The form for 'External Energy 1' includes the following fields:


- Name: IKA Cotton Thread
- Gross Calorific Value [J/g]: (radio button selected, circled in red with '1')
- External Energy (fix) [J]: 50 (radio button selected, circled in red with '2')
- Note: Auto generated entry (circled in red with '3')

Buttons on the right include: New External Energy, Delete External Energy, Change External Energy, and Reset External Energy.

External Energy Name	Gross Heat Value	External Energy
IKA Cotton Thread	0	50

Buttons at the bottom: Print, OK

3. Select your kind of external energy:
 1. Gross Calorific Value: an additional indication of the weight is necessary when setting up a new measurement
 2. External Energy (fix): predefined fix energy value
 3. It is possible to add a note to the external energy that's being created, e.g., creation date



4. Confirm by clicking “Add External Energy”

Settings

External Energy 1 External Energy 2 External Energy 3

User Management

External Energies

Calorimeter

Vessels

C-Values

System Settings

Groupings

Name:

Gross Calorific Value [J/g]:

External Energy (fix) [J]:

Note:

New External Energy

Delete External Energy

Change External Energy

Reset External Energy

External Energy Name	Gross Heat Value	External Energy
Ext. Energy 1	22450	0
IKA Cotton Thread	0	50

Print OK

5. Now it is possible to select this external energy as "Combustion Aid 1" in the new measurement screen and to add a weight.

New measurement

Messung

Measurement name:

User:

Calorimeter:

Vessel:

Weight:

Group:

Calibration

Simulation

User Defined Reference Combustion Value

Demo Mode

Combustion Aid 1

Combustion Aid 2

Combustion Aid 3


Type:	Hot:	Weight:	QExt:
Ext. Energy 1	22450	0,5	11225

Sum External Energies:


Information:

Clear OK Cancel

5.3 C-Values



1. Select a calorimeter
2. Select a vessel
3. Select a measurement method
4. Press the "Start Searching" button



5. The found calibrations are now listed.
6. Select the calibrations for calculating the c-value. The selected calibrations are shown in the graph.

The screenshot shows the 'Settings' window with the 'C-Values' section selected in the left sidebar. A table lists two probes with their respective C-values. To the right, a summary of statistics is displayed, including an average C-value of 5186.3. A red circle highlights the checkboxes for both probes in the table. A red arrow points from the 'Average' value in the summary to the 'C-Value' field in the 'Vessel Informations' window shown in the next screenshot. Below the table is a line graph showing the C-value over time for the two probes, with a mean line and upper/lower sigma limits.

Probe Name	Start Date	Result C-Value
1211291116	11/29/2012 11:16:24 ...	4947.3
1211291122	11/29/2012 11:22:24 ...	5425.3

Summary Statistics:
 Average: 5186.3
 RSD: 100.2121
 Calculated σ [J]: 0
 Current σ [%]: 0
 Change σ

7. Click „Save Values“ to set the calculated average as new c-value.
8. The new c-value can also be found in the vessels settings form.

The screenshot shows the 'Settings' window with the 'Vessels' section selected in the left sidebar. The 'Vessel Informations' form is displayed, showing various fields for vessel configuration. The 'C-Value' field is set to 5186, which is circled in red. A red arrow points from the 'Average' value in the first screenshot to this 'C-Value' field. Below the form is a table listing vessel information, with the 'C-Value' column also circled in red.

Vessel Informations:

Vessel-ID/RFID: 1
 Name: C1 Vessel No 1
 Optical ID:
 Calorimeter: C1 No 1 (Isoperibol 22°C | COM3)
 Working Mode: Isoperibol 22°C
 C-Value: 5186
 Ignitions: 12
 Note:

Vessel Name	Measurement Mode	Vessel-ID/RFID	Optical ID	C-Value	Ignition Counter
C1 Vessel No 1	Isoperibol 22°C	1		5186	12

5.3.1 Comparing C-Values

It is possible to check if the c-values are inside a tolerance limit (σ - sigma).


The screenshot shows the 'Settings' window with the following details:

- Configuration:** Calorimeter: C1 No 1 (Isoperibol 22°C | COM3), Vessel: C1 Vessel No 1 (ID: 1), Measurement Method: Isoperibol 22°C. A 'Start Searching' button is present.
- C-Values Table:**

	Probe Name	Start Date	Result C-Value
<input checked="" type="checkbox"/>	1211291116	11/29/2012 11:16:24 ...	4947.3
<input checked="" type="checkbox"/>	1211291122	11/29/2012 11:22:24 ...	5425.3
- Statistics:** Average: 5186.3, RSD: 100.2121, Calculated σ [J]: 259.315, Current σ [%]: 5.
- Controls:** A 'Change σ ' checkbox (1) is checked, with an input field (2) containing '5'. A 'Save Values' button (3) is located below.
- Graph:** A line graph (4) plots C-Value (y-axis, 4800-5600) against Time (x-axis, 29-Nov-12). It shows a blue line connecting two data points, with horizontal lines for 'Upper σ ', 'Mean', and 'Lower σ '.

1. Select the “Change σ ” checkbox
2. Enter a value into the appearing input field.
3. Press the “Save Values” button to apply the input value
4. Check the upper and lower limits in the graph.

5.4 System Settings




1. Select the evaluation norm to calculate the result of a measurement.
2. The calculation unit.
3. The result value used for the calculations in the “Analysis” tab in the main screen.
4. The application language.
5. The Reference Combustion Value of Benzoic Acid.
6. Possible to connect a balance directly with the PC. If a port is selected, the port setting input fields get activated and the port test button.
7. It is possible to exclude unneeded COM-Ports. If a port cannot be opened or is not reachable at startup it will be added to this list.
If you want to re-activate a port, remove it from the input field and it will be used on next application start.
8. Reset the user settings to initial values.
9. If you want to save changes of any input field(s), press the “Save Changes” button to take over the new value(s).

5.4.1 Grid Columns 1 – 3

Set the visible columns at the main screen. If you want to hide or add a column, this is possible in the 3 “Grid Columns” tabs.

To take over the changes on a tab, press on the “Save Changes” button.



5.4.2 Database Settings

On this form you can change the access to the underlying database. This tab is only visible if you selected the SQL Data Storage method.

WARNING: change these values only if you know exactly what you are doing!!!


The screenshot shows the 'Database Settings' window with the following details:

- Server Name:** SL-11-030
- Database Name:** CalDBTest
- Server Login:**
- User Name:** [Empty text box]
- Password:** [Empty text box]
- Confirm password:** [Empty text box]
- Error Message:** Passwords not equal!
- URL:** <http://www.connectionstrings.com/sql-server-2008>
- Buttons:** Get local servers, Get network servers, Request Databases, Check Database Connection, Save Changes, OK


5.5 Groupings

You can create various groups, to have an additional way to sort your measurements.


E.g.: If you have samples which should be tested on different days, you can add them to a group so that they can be found, displayed and analyzed much easier on the main screen.




1. Press "New Group".
2. Enter Name.




3. Now you can add measurements from the main screen to this group.
4. Open measurement on the main screen by double clicking.




5. Measurement opens and Group can be selected



6. Select Group and press on the "Close" button



7. The selected group name is shown in the Details tab and in the Group Name column, if it is activated in the Settings.



8. Now you can use the group name as a filtering option when searching for measurements. Just select the group you want to display and click on "Update Measurements".

6 Detail tabs

The screenshot shows the IKA CalWin software interface. At the top, there is a menu bar with options: New Measurement, Settings, Export, Print, and Help. The main window is divided into several sections:

- Left Panel:** Contains navigation menus for Period, Calorimeters, User, and Groups. The 'Imported C1' group is selected.
- Selected Measurements Table:** A table with columns: Calorimeter Name, Probe Name, Vessel Name, Weight [g], Result C-Value, Result (Ho), Vessel C-Value, and User Name. Two rows are visible, both for 'C1 No 1'.
- Device View:** Shows 'C1 No 1' with a 'Waiting...' status and a 'Start' button.
- Details Panel (circled in red):** Provides a detailed view of the selected measurement, including:
 - Ext. Energy 3 Calculated Value: --
 - Ignition Energy: 9
 - Delta T: 5.8913
 - Vessel C-Value: 4000
 - State: completed
 - Group Name: Imported C1
 - Measurement Duration: 00:00:00
 - Creation date: 11/28/2012 12:15:30 PM
 - Ref. Comb. Value: 26457
 - Calorimeter Type: C1
 - Measurement Mode: Isoperibol 22°C

6.1 Details


Get an overview of the selected measurement.

This screenshot shows a close-up of the 'Details' tab in the IKA CalWin software. It displays the following information:

- Probe Name:** 1211281215
- Weight [g]:** 1.22
- Result (Ho):** 19268
- User Name:** notSet
- Calorimeter Name:** C1 No 1
- Vessel Name:**
- Start Date:** 11/28/2012 12:15:30 PM
- External Energies:**
 - Ext. Energy 1 Calculated Value: 50
 - Ext. Energy 2 Calculated Value: --
 - Ext. Energy 3 Calculated Value: --
 - Ignition Energy: 9

6.2 Graph


See the temperature curve of a selected measurement. This doesn't work for imported measurements!



6.3 Analysis

Select measurements using the check box and compare them in the "Analysis" tab.


Selected Measurements:											
<input type="checkbox"/>		Calorimeter Name	Probe Name	Vessel Name	Weight [g]	Result C-Value	Result (Ho)	Vessel C-Value	User Name	Delta T	Ext. Er
<input checked="" type="checkbox"/>	✓	C1 No 1	1211281209		1.11	0	21108	4000	notSet	5.8721	50
<input checked="" type="checkbox"/>	✓	C1 No 1	1211281215		1.22	0	19268	4000	notSet	5.8913	50
<input type="checkbox"/>	✓	C1 No 1	1211291116		1.1	4947.3	0	4000	notSet	5.8943	50
<input type="checkbox"/>	✓	C1 No 1	1211291122		1.2	5425.3	0	4000	notSet	5.8627	50




NOTE: these values are just an example for visual representation of σ .

6.4 Notes

It's possible to add some notes to a selected measurement



1. Enter note to add
2. Click "Add Note". The Note is now added to the Note history.



INFORMATION: Notes can't be removed!

6.5 Printing/Export Details

Every tab can be printed or exported to EXCEL. Just click on the button on the top right corner. The print button opens a printer dialog and the export button creates an EXCEL file under: C:\user documents \CaWin\Export that contains the exported data.

