

## Tip compatibility chart

Pipette brand	Pipette model	Tip 10 $\mu$ l bag Ident. No. 0020114634	Tip 200 $\mu$ l bag Ident. No 0020114635	Tip 1000 $\mu$ l bag Ident. No 0020114636
IKA	<b>PETTE fix</b> single-channel, fixed volume			
	5 $\mu$ l	X		
	10 $\mu$ l	X		
	20 $\mu$ l		X	
	25 $\mu$ l		X	
	50 $\mu$ l		X	
	100 $\mu$ l		X	
	200 $\mu$ l		X	
	250 $\mu$ l			X
	500 $\mu$ l			X
	1000 $\mu$ l			X
	<b>PETTE vario</b> single-channel, variable volume			
	0.1 – 2 $\mu$ l	X		
	0.5 – 10 $\mu$ l	X		
	2 – 20 $\mu$ l		X	
	10 – 100 $\mu$ l		X	
	20 – 200 $\mu$ l		X	
	100 – 1000 $\mu$ l			X
	0.5 – 5 ml			
	1 – 10 ml			
	<b>PETTE multi</b> 8-channel, variable volume			
	8 x 0.5 - 10 $\mu$ l	X		
	8 x 10 - 100 $\mu$ l		X	
	8 x 20 - 200 $\mu$ l		X	
	<b>PETTE multi</b> 12-channel, variable volume			
	12 x 0.5 - 10 $\mu$ l	X		
	12 x 10 - 100 $\mu$ l		X	
	12 x 20 - 200 $\mu$ l		X	

## Tip compatibility chart

Pipette brand	Pipette model	Tip 10 ul bag Ident. No. 0020114634	Tip 200 ul bag Ident. No 0020114635	Tip 1000 ul bag Ident. No 0020114636
<b>EPPENDORF</b>				
	<b>Reference</b>			
	0.5 - 10 µl	x		
	10 - 100 µl		x	
	50 - 200 µl		x	
	200 - 1000 µl			x
	<b>Research</b>			
	0.5 - 10 µl	x		
	2 - 20 µl		x	
	10 - 100 µl		x	
	20 - 200 µl		x	
	100 - 1000 µl			x
8-channel	0.5 - 10 µl		x	
	<b>Reference 2</b> single-channel, variable volume			
	0.5 - 10 µl	x		
	2 - 20 µl		x	
	10 - 100 µl		x	
	20 - 200 µl		x	
	100 - 1000 µl			x
	0.25 - 2.5 ml			
	0.5 - 5 ml			
	1 - 10 ml			
8-channel	10 - 100 µl		x	
	<b>Research Plus</b> single-channel, variable volume			
	0.5 - 10 µl	x		
	2 - 20 µl	x		
	10 - 100 µl		x	
	20 - 200 µl		x	
	100 - 1000 µl			x
	0.5 - 5 ml			
	1 - 10 ml			
8-channel	0.5 - 10 µl	x		
8-channel	10 - 100 µl		x	
12-channel	0.5 - 10 µl	x		
12-channel	10 - 100 µl		x	

Please note that the above compatibility table is for guidance only.  
 IKA is not liable for technical modifications or changes to the listed pipettes.  
 The specifications (= precision and accuracy) have not been tested.

## Tip compatibility chart

Pipette brand	Pipette model	Tip 10 ul bag Ident. No. 0020114634	Tip 200 ul bag Ident. No 0020114635	Tip 1000 ul bag Ident. No 0020114636	
<b>SARTORIUS</b>	<b>Proline</b>				
		0.5 - 10 µl	X		
		10 - 100 µl		X	
		20 - 200 µl		X	
		100 - 1000 µl			X
		1000 - 5000 µl			
	<b>mLine</b>				
		0.5 - 10 µl	X		
		10 - 100 µl		X	
		20 - 200 µl		X	
		100 - 1000 µl			X
		500 - 5000 µl			
	<b>GILSON</b>	<b>Pipetman Classic / G / L - traditional shaft</b>			
		P10	1 - 10 µl	X	
P20		2 - 20 µl	X		
P100		10 - 100 µl		X	
P200		50 - 200 µl		X	
P1000		100 - 1000 µl			X
P5000		1 - 5 ml			
P10ML		1 - 10 ml			
<b>Pipetman Ultra</b>					
U10		1 - 10 µl	∅*		
U100		10 - 100 µl		X	
U200		20 - 200 µl		X	
U1000		200 - 1000 µl			X
U5000		1 - 5 ml			
U10000		1 - 10 ml			
<b>Pipetman Neo</b>					
		1 - 10 µl	X		
		2 - 20 µl	∅*		
		10 - 100 µl		X	
	20 - 200 µl		X		
	100 - 1000 µl			X	

Please note that the above compatibility table is for guidance only. IKA is not liable for technical modifications or changes to the listed pipettes. The specifications (= precision and accuracy) have not been tested.

\* Tip ejection is not possible. It has to be manually removed.

## Tip compatibility chart

Pipette brand	Pipette model	Tip 10 ul bag Ident. No. 0020114634	Tip 200 ul bag Ident. No 0020114635	Tip 1000 ul bag Ident. No 0020114636
SOCOREX	<b>Acura 825</b> single channel digital			
	0.5 - 10 µl	x		
	2 - 20 µl			
	10 - 100 µl		x	
	20 - 200 µl		x	
	100 - 1000 µl			
	<b>Acura manual 835</b>			
	0.2 - 2 ml			
	0.5 - 5 ml			
	1 - 10 ml			
	<b>Acura manual XS</b>			
	0.5 - 10 µl	x		
	2 - 20 µl			
	10 - 100 µl		x	
	20 - 200 µl		x	
	100 - 1000 µl			x
	<b>Calibra Digital 822</b>			
	1 - 10 µl	x		
	2 - 20 µl			
	10 - 100 µl			
	20 - 200 µl		x	
	100 - 1000 µl			x
	<b>Calibra Digital 832</b>			
	0.2 - 2 ml			
	1 - 10 ml			

Please note that the above compatibility table is for guidance only.  
 IKA is not liable for technical modifications or changes to the listed pipettes.  
 The specifications (= precision and accuracy) have not been tested.

## Tip compatibility chart

Pipette brand	Pipette model	Tip 10 ul bag Ident. No. 0020114634	Tip 200 ul bag Ident. No 0020114635	Tip 1000 ul bag Ident. No 0020114636
<b>BRAND</b>	<b>Transferpette S</b>			
	0.5 - 10 µl	X		
	2 - 20 µl		X	
	10 - 100 µl		X	
	20 - 200 µl		X	
	100 - 1000 µl			X
	0.5 - 5 ml			
	1 - 10 ml			
	<b>Transferpette electronic</b>			
	0.5 - 10 µl	X		
	2 - 20 µl		X	
	10 - 100 µl		X	
	20 - 200 µl		X	
	100 - 1000 µl			X
200 - 5000 µl				
<b>THERMO FINNPIPETTE</b>	<b>Finnpipette digital pipettes and F1 / F2 / F3</b>			
	0.2 - 2 µl	X		
	0.5 - 10 µl	X		
	10 - 100 µl		X	
	20 - 200 µl		X	
	200 - 1000 µl			X
	1 - 5 ml			
2 - 10 ml				

Please note that the above compatibility table is for guidance only.  
 IKA is not liable for technical modifications or changes to the listed pipettes.  
 The specifications (= precision and accuracy) have not been tested.