



designed for scientists

Maximum power with minimum footprint

/// The smart mini disperser series

NEW

T 18 mini digital
ULTRA-TURRAX®
Ident. No. 0025005828

T 25 mini control
ULTRA-TURRAX®
Ident. No. 0025006650

Meet the new mini dispersers

Introducing the latest in high-shear dispersing technology: the **T 18 mini digital** and **T 25 mini control**. These new additions to the IKA ULTRA-TURRAX® family are designed for precision dispersing, emulsifying, and homogenizing, offering high performance in a compact, innovative way.

T 18 mini digital & T 25 mini control

Redefine your lab's efficiency and precision with IKA's new mini homogenizers. Compact, powerful, and ready to meet your most demanding requirements.

ADVANTAGES

T 18 mini digital / T 25 mini control

- › **Compact and Space-Saving Design:** up to 40% smaller footprint than comparable models, fits into even the smallest lab spaces.
- › **Cleanroom Ready:** Brushless DC motor ensures dust-free and quiet operation.
- › **Clear TFT Display:** Easy-to-read screen for speed and timer settings.
- › **Smart Connectivity:** Full range of interfaces enables software updates, communication with other instruments and integration with labworldsoft® software for automation and data logging.
- › **Unattended Operation:** Timer and counter functions ensures safe working and reproducible test results.
- › **Perfect for your low-viscosity homogenization tasks:** With its motor power of 160 W, the devices are ideal for any simple homogenizing task up to 1.5 liters and support viscosities up to 2000 mPas.

ADDITIONAL ADVANTAGES

T 25 mini control

- › **Additional Safety for Samples:** The patented temperature sensing function in the dispersing tools monitors medium temperature in real-time. Safety temperature settings protect temperature-sensitive samples from degradation.
- › **Easy to clean and time-saving:** Compatible with "easy-clean" dispersing tools for quick and simple cleaning.

up to 40%
smaller



For more information, visit www.ika.com or contact your local IKA Sales Team.

IKA-Werke GmbH & Co. KG

Janke & Kunkel-Straße 10, 79219 Staufen, Germany

Phone: +49 7633 831-0, eMail: sales@ika.de



www.ika.com



IKAworlwide // #lookattheblue