

LAB ESSENTIALS FOR OIL & GAS INDUSTRY



Magnetic Stirrers

IKA Plate (RCT digital)

> General lab mixing and heating



Overhead Stirrer

EUROSTAR 60 control

> Large volume mixing, blending and solutions



Disperser | Homogenizer

T 25 easy clean control ULTRA-TURRAX®

> High-shear emulsification and mixing processes



Viscometers

ROTAVISC me-vi Complete

> Determine the viscosity of liquids in all areas of application ranging from laboratory to quality control



Calorimeters

C 6000 Oxygen Combustion Calorimeter

> Used in quality control, determination of degree of efficiency of engines, and heavy fuel oil testing (HFO), etc.
> Complies to the standard test methods for heat of combustion of ASTM D 240 and ASTM D 4809 (precision method).
> Complies to DIN 51900 on solid and liquid fuels



Pipettes & Accessories

IKA PETTE fix and vario

> Precise dosing, measuring and preparation of liquid samples

