

LAB ESSENTIALS FOR PAINTING AND COATING INDUSTRY

ROTAISC SBS me-vi Complete

Viscosity analysis, quality control, R&D

- > High robustness, high sample throughput, precise measurement, operation in harsher environments.



EUROSTAR 100 digital

Large volume mixing, blending and solutions

- > Highly viscous applications and intensive mixing with automatic cut-off under anti-stall or overload conditions, constant speed throughout changing viscosities.



C-MAG HS 7 control

Mixing and heating with precise temperature control

- > Ceramic heating plate with excellent chemical resistance, stirring quantities of up to 10 l, heating temperature range 50 °C – 550 °C.



ROLLER 10 digital

Gentle motion and uniformity for solid and liquid suspensions and viscous samples

- > Smooth rocking and rolling action at infinitely adjustable speeds, compatible with all conventional tubes and cylindrical bottles.



MultiDrive basic

Coarse & fine crushing and milling

- > Crushing performance guaranteed by a combination of variable rotational speeds, interval operation, integrated cooling and a variety of vessels.



IKA OVEN 125 basic dry – glass

Tempering, dry-aging

- > Fast heat up times and a precise temperature control guarantee reproducible test results.

