

Tip compatibility chart

Pipette brand	Pipette model	IKA Tip 10 µl bag Ident. No. 0020114634	IKA Tip 200 µl bag Ident. No. 0020114635	IKA Tip 1000 µl bag Ident. No. 0020114636
		IKA Tip 10 µl box Ident. No. 0020116857	IKA Tip 200 µl box Ident. No. 0020116858	IKA Tip 1000 µl box Ident. No. 0020116859
IKA	PETTE fix single-channel, fixed volume			
	5 µl	x		
	10 µl	x		
	20 µl		x	
	25 µl		x	
	50 µl		x	
	100 µl		x	
	200 µl		x	
	250 µl			x
	500 µl			x
	1000 µl			x
	PETTE vario single-channel, variable volume			
	0.1 – 2 µl	x		
	0.5 – 10 µl	x		
	2 – 20 µl		x	
	10 – 100 µl		x	
	20 – 200 µl		x	
	100 – 1000 µl			x
	0.5 – 5 ml			
	1 – 10 ml			
	PETTE multi 8-channel, variable volume			
	8 x 0.5 - 10 µl	x		
	8 x 10 - 100 µl			
	8 x 20 - 200 µl			
	PETTE multi 12-channel, variable volume			
	12 x 0.5 - 10 µl	x		
	12 x 10 - 100 µl			
	12 x 20 - 200 µl			
	Electronic Pipette single-channel, variable volume			
	0.5 - 10 µl	x		
	5 - 100 µl		x	
	10 - 200 µl		x	

Tip compatibility chart

Pipette brand	Pipette model	IKA Tip 10 µl bag Ident. No. 0020114634	IKA Tip 200 µl bag Ident. No. 0020114635	IKA Tip 1000 µl bag Ident. No. 0020114636
		IKA Tip 10 µl box Ident. No. 0020116857	IKA Tip 200 µl box Ident. No. 0020116858	IKA Tip 1000 µl box Ident. No. 0020116859
EPPENDORF	Reference			
	0.5 - 10 µl	x		
	10 - 100 µl		x	
	50 - 200 µl		x	
	200 - 1000 µl			x
	Research			
	0.5 - 10 µl	x		
	2 - 20 µl		x	
	10 - 100 µl		x	
	20 - 200 µl		x	
	100 - 1000 µl			x
	8-channel 0.5 - 10 µl		x	
	Reference 2 single-channel, variable volume			
	0.5 - 10 µl	x		
	2 - 20 µl		x	
	10 - 100 µl		x	
	20 - 200 µl		x	
	100 - 1000 µl			x
	0.25 - 2.5 ml			
	0.5 - 5 ml			
	1 - 10 ml			
	8-channel 10 - 100 µl		x	
	Research Plus single-channel, variable volume			
	0.5 - 10 µl	x		
	2 - 20 µl	x		
	10 - 100 µl		x	
	20 - 200 µl		x	
	100 - 1000 µl			x
	0.5 - 5 ml			
	1 - 10 ml			
	8-channel 0.5 - 10 µl	x		
	8-channel 10 - 100 µl		x	
	12-channel 0.5 - 10 µl	x		
	12-channel 10 - 100 µl		x	

Please note that the above compatibility table is for guidance only. IKA is not liable for technical modifications or changes to the listed pipettes. The specifications (= precision and accuracy) have not been tested.



designed for scientists

Tip compatibility chart

Pipette brand	Pipette model	IKA Tip 10 µl bag Ident. No. 0020114634	IKA Tip 200 µl bag Ident. No. 0020114635	IKA Tip 1000 µl bag Ident. No. 0020114636
		IKA Tip 10 µl box Ident. No. 0020116857	IKA Tip 200 µl box Ident. No. 0020116858	IKA Tip 1000 µl box Ident. No. 0020116859
SARTORIUS	Proline			
	0.5 - 10 µl	x		
	10 - 100 µl		x	
	20 - 200 µl		x	
	100 - 1000 µl			x
	1000 - 5000 µl			
	mLine			
	0.5 - 10 µl	x		
	10 - 100 µl		x	
	20 - 200 µl		x	
	100 - 1000 µl			x
	500 - 5000 µl			
GILSON	Pipetman Classic / G / L - traditional shaft			
	P10	1 - 10 µl	x	
	P20	2 - 20 µl	x	
	P100	10 - 100 µl	x	
	P200	50 - 200 µl	x	
	P1000	100 - 1000 µl		x
	P5000	1 - 5 ml		
	P10ML	1 - 10 ml		
	Pipetman Ultra			
	U10	1 - 10 µl	ø*	
	U100	10 - 100 µl	x	
	U200	20 - 200 µl	x	
	U1000	200 - 1000 µl		x
	U5000	1 - 5 ml		
	U10000	1 - 10 ml		
	Pipetman Neo			
		1 - 10 µl	x	
		2 - 20 µl	ø*	
		10 - 100 µl	x	
		20 - 200 µl	x	
		100 - 1000 µl		x

* Tip ejection is not possible. It has to be manually removed.



designed for scientists

Tip compatibility chart

Pipette brand	Pipette model	IKA Tip 10 µl bag Ident. No. 0020114634	IKA Tip 200 µl bag Ident. No. 0020114635	IKA Tip 1000 µl bag Ident. No. 0020114636
		IKA Tip 10 µl box Ident. No. 0020116857	IKA Tip 200 µl box Ident. No. 0020116858	IKA Tip 1000 µl box Ident. No. 0020116859
SOCOREX	Acura 825 single channel digital			
	0.5 - 10 µl	x		
	2 - 20 µl			
	10 - 100 µl		x	
	20 - 200 µl		x	
	100 - 1000 µl			
	Acura manual 835			
	0.2 - 2 ml			
	0.5 - 5 ml			
	1 - 10 ml			
	Acura manual XS			
	0.5 - 10 µl	x		
	2 - 20 µl			
	10 - 100 µl		x	
	20 - 200 µl		x	
	100 - 1000 µl			x
	Calibra Digital 822			
	1 - 10 µl	x		
	2 - 20 µl			
	10 - 100 µl			
	20 - 200 µl		x	
	100 - 1000 µl			x
	Calibra Digital 832			
	0.2 - 2 ml			
	1 - 10 ml			

Please note that the above compatibility table is for guidance only. IKA is not liable for technical modifications or changes to the listed pipettes. The specifications (= precision and accuracy) have not been tested.

Tip compatibility chart

Pipette brand	Pipette model	IKA Tip 10 µl bag Ident. No. 0020114634	IKA Tip 200 µl bag Ident. No. 0020114635	IKA Tip 1000 µl bag Ident. No. 0020114636
		IKA Tip 10 µl box Ident. No. 0020116857	IKA Tip 200 µl box Ident. No. 0020116858	IKA Tip 1000 µl box Ident. No. 0020116859
BRAND	Transferpette S			
	0.5 - 10 µl	x		
	2 - 20 µl		x	
	10 - 100 µl		x	
	20 - 200 µl		x	
	100 - 1000 µl			x
	0.5 - 5 ml			
	1 - 10 ml			
	Transferpette electronic			
	0.5 - 10 µl	x		
	2 - 20 µl		x	
	10 - 100 µl		x	
	20 - 200 µl		x	
	100 - 1000 µl			x
	200 - 5000 µl			
	Finnpipette digital pipettes and F1 / F2 / F3			
	0.2 - 2 µl	x		
	0.5 - 10 µl	x		
	10 - 100 µl		x	
	20 - 200 µl		x	
	200 - 1000 µl			x
	1 - 5 ml			
	2 - 10 ml			
THERMO FINNPIPETTE				

Please note that the above compatibility table is for guidance only. IKA is not liable for technical modifications or changes to the listed pipettes. The specifications (= precision and accuracy) have not been tested.