

Tip compatibility chart

Pipette brand	Pipette model	IKA Tip 10 µl bag Ident. No. 0020114634	IKA Tip 200 µl bag Ident. No. 0020114635	IKA Tip 1000 µl bag Ident. No. 0020114636
		IKA Tip 10 µl box Ident. No. 0020116857	IKA Tip 200 µl box Ident. No. 0020116858	IKA Tip 1000 µl box Ident. No. 0020142071
IKA	PETTE fix single-channel, fixed volume			
	5 µl	x		
	10 µl	x		
	20 µl		x	
	25 µl		x	
	50 µl		x	
	100 µl		x	
	200 µl		x	
	250 µl			x
	500 µl			x
	1000 µl			x
	PETTE vario single-channel, variable volume			
	0.1 – 2 µl	x		
	0.5 – 10 µl	x		
	2 – 20 µl		x	
	10 – 100 µl		x	
	20 – 200 µl		x	
	100 – 1000 µl			x
	0.5 – 5 ml			
	1 – 10 ml			
	PETTE multi 8-channel, variable volume			
	8 x 0.5 – 10 µl	x		
	8 x 10 – 100 µl		x	
	8 x 20 – 200 µl		x	
	PETTE multi 12-channel, variable volume			
	12 x 0.5 – 10 µl	x		
	12 x 10 – 100 µl		x	
	12 x 20 – 200 µl		x	

Tip compatibility chart

Pipette brand	Pipette model	IKA Tip 10 µl bag Ident. No. 0020114634	IKA Tip 200 µl bag Ident. No. 0020114635	IKA Tip 1000 µl bag Ident. No. 0020114636
		IKA Tip 10 µl box Ident. No. 0020116857	IKA Tip 200 µl box Ident. No. 0020116858	IKA Tip 1000 µl box Ident. No. 0020142071
EPPENDORF	Reference			
	0.5 – 10 µl	x		
	10 – 100 µl		x	
	50 – 200 µl		x	
	200 – 1000 µl			x
	Research			
	0.5 – 10 µl	x		
	2 – 20 µl		x	
	10 – 100 µl		x	
	20 – 200 µl		x	
	100 – 1000 µl			x
	8-channel 0.5 – 10 µl		x	
	Reference 2 single-channel, variable volume			
	0.5 – 10 µl	x		
	2 – 20 µl		x	
	10 – 100 µl		x	
	20 – 200 µl		x	
	100 – 1000 µl			x
	0.25 – 2.5 ml			
	0.5 – 5 ml			
	1 – 10 ml			
	8-channel 10 – 100 µl		x	
	Research Plus single-channel, variable volume			
	0.5 – 10 µl	x		
	2 – 20 µl	x		
	10 – 100 µl		x	
	20 – 200 µl		x	
	100 – 1000 µl			x
	0.5 – 5 ml			
	1 – 10 ml			
	8-channel 0.5 – 10 µl	x		
	8-channel 10 – 100 µl		x	
	12-channel 0.5 – 10 µl	x		
	12-channel 10 – 100 µl		x	

Tip compatibility chart

Pipette brand	Pipette model	IKA Tip 10 µl bag Ident. No. 0020114634	IKA Tip 200 µl bag Ident. No. 0020114635	IKA Tip 1000 µl bag Ident. No. 0020114636
		IKA Tip 10 µl box Ident. No. 0020116857	IKA Tip 200 µl box Ident. No. 0020116858	IKA Tip 1000 µl box Ident. No. 0020142071
SARTORIUS	Proline			
	0.5 – 10 µl	x		
	10 – 100 µl		x	
	20 – 200 µl		x	
	100 – 1000 µl			x
	1000 – 5000 µl			
	mLine			
	0.5 – 10 µl	x		
	10 – 100 µl		x	
	20 – 200 µl		x	
	100 – 1000 µl			x
	500 – 5000 µl			
GILSON	Pipetman Classic / G / L - traditional shaft			
	P10 1 – 10 µl	x		
	P20 2 – 20 µl	x		
	P100 10 – 100 µl		x	
	P200 50 – 200 µl		x	
	P1000 100 – 1000 µl			x
	P5000 1 – 5 ml			
	P10ML 1 – 10 ml			
	Pipetman Ultra			
	U10 1 – 10 µl	Ø*		
	U100 10 – 100 µl		x	
	U200 20 – 200 µl		x	
	U1000 200 – 1000 µl			x
	U5000 1 – 5 ml			
	U10000 1 – 10 ml			

* Tip ejection is not possible.
It has to be manually removed.



designed for scientists

Tip compatibility chart

Pipette brand	Pipette model	IKA Tip 10 µl bag Ident. No. 0020114634	IKA Tip 200 µl bag Ident. No. 0020114635	IKA Tip 1000 µl bag Ident. No. 0020114636
		IKA Tip 10 µl box Ident. No. 0020116857	IKA Tip 200 µl box Ident. No. 0020116858	IKA Tip 1000 µl box Ident. No. 0020142071
GILSON	Pipetman Neo			
	1 – 10 µl	x		
	2 – 20 µl	Ø*		
	10 – 100 µl		x	
	20 – 200 µl		x	
	100 – 1000 µl			x
SOCOREX	Acura 825 single channel digital			
	0.5 – 10 µl	x		
	2 – 20 µl			
	10 – 100 µl		x	
	20 – 200 µl		x	
	100 – 1000 µl			
	Acura manual 835			
	0.2 – 2 ml			
	0.5 – 5 ml			
	1 – 10 ml			
	Acura manual XS			
	0.5 – 10 µl	x		
	2 – 20 µl			
	10 – 100 µl		x	
	20 – 200 µl		x	
	100 – 1000 µl			x
	Calibra Digital 822			
	1 – 10 µl	x		
	2 – 20 µl			
	10 – 100 µl			
	20 – 200 µl		x	
	100 – 1000 µl			x
	Calibra Digital 832			
	0.2 – 2 ml			
	1 – 10 ml			

Please note that the above compatibility table is for guidance only. IKA is not liable for technical modifications or changes to the listed pipettes. The specifications (= precision and accuracy) have not been tested.

Tip compatibility chart

Pipette brand	Pipette model	IKA Tip 10 µl bag Ident. No. 0020114634	IKA Tip 200 µl bag Ident. No. 0020114635	IKA Tip 1000 µl bag Ident. No. 0020114636
		IKA Tip 10 µl box Ident. No. 0020116857	IKA Tip 200 µl box Ident. No. 0020116858	IKA Tip 1000 µl box Ident. No. 0020142071
BRAND	Transferpette S			
	0.5 – 10 µl	x		
	2 – 20 µl		x	
	10 – 100 µl		x	
	20 – 200 µl		x	
	100 – 1000 µl			x
	0.5 – 5 ml			
	1 – 10 ml			
	Transferpette electronic			
	0.5 – 10 µl	x		
	2 – 20 µl		x	
	10 – 100 µl		x	
	20 – 200 µl		x	
	100 – 1000 µl			x
	200 – 5000 µl			
	Finnpipette digital pipettes and F1 / F2 / F3			
	0.2 – 2 µl	x		
	0.5 – 10 µl	x		
	10 – 100 µl		x	
	20 – 200 µl		x	
	200 – 1000 µl			x
	1 – 5 ml			
	2 – 10 ml			
THERMO FINNPIPETTE				

Please note that the above compatibility table is for guidance only. IKA is not liable for technical modifications or changes to the listed pipettes. The specifications (= precision and accuracy) have not been tested.